


Application #1

Single Bench Fume Hood



In this Sequence:

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #1: Single Bench Fume Hood

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LEDs capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood and perpendicular to the sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an equal

as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If velocity continues to fall below setpoint, the controller generates a local alarm (as described in the alarming section) after an adjustable time delay. When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If velocity continues to rise above setpoint, the controller generates a local alarm after an adjustable time delay.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Such devices usually consist of emergency panic pushbuttons located at the exits of a space. Initiation of either device positions the hood exhaust damper to an user-defined fixed position. In this mode the controller display indicates that the hood is in emergency mode while the audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms both locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|---------|-----------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |


percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned

Sequence of Operation



In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms can be silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that mimics the operation of the built-in audible alarm.

3



NOTE:
SEE DRAWING A950051
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.


TYPICAL FUME HOOD LAB
APPLICATION H1

| | | | | |
|--------------------------------------|-----------|------------|---------|-------|
| DRAWN LIP | ENG'D. BY | CHECKED BY | DATE | SHEET |
| | | | 9-12-94 | 1 |
| TYPICAL FUME HOOD LAB APPLICATION H1 | | | | |
| DRAWING NO. A940041 | | | | |

BENCH TYPE FUME HOOD

Application #2

Single Walk-in or Distillation
Fume Hood



In this Sequence:

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #2: Single Walk-in or Distillation Fume Hood

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening with a space pressure primary sensor. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading acts as the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or

the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller with a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop 20 times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If velocity continues to rise above setpoint, the controller generates a local alarm after an adjustable time delay.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |


remotely. The controller's PID output is calculated by

Sequence of Operation



In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module.. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.

7



NOTE:
SEE DRAWING A950052
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.


| | | | |
|---------------|--|------------------------------|--|
| DRAWN LIP | | AUTO-FLOW | |
| ENGR. BY | | BY TECHNOLOGY DEPT. 12/19/82 | |
| CHECKED BY | | SINGLE WALK-IN HOOD | |
| DATE: 9-13-95 | | FLOW DIAGRAM | |
| SHEET | | DRAWING NO. A910042 | |

TYPICAL FUME HOOD LAB
APPLICATION #2

WALK-IN OR DISTILLATION
TYPE FUME HOOD

Application #3

Single Bench Fume Hood with
Air Flow Tracking



In this Sequence:

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #3: Single Bench Fume Hood with Air Flow Tracking

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LEDs capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face

velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for use by the laboratory pressure controller in calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

Sequence of Operation




| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period.

After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.

[illegible]




NOTE:
SEE DRAWING A950053
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

BENCH TYPE FUME HOOD
W/EXHAUST FLOW MEASUREMENT

| | | | | | |
|---------------------------------------|--|---------------|--|----------------------------------------|--|
| TYPICAL FUME HOOD LAB APPLICATION # 3 | | DRAWN LIP | | AUTO-FLOW DO NOT SCALE, SEE P. 1592 | |
| | | ENGR: RVF | | SINGLE BENCH | |
| | | CHKD: RVF | | AIR FLOW TRACKING FLOW DIAGRAM | |
| | | DATE: 8-13-99 | | | |
| | | SHEET 3 | | DRAWING NO. A940043 | |

Application #4

Single Walk-in or Distillation
Fume Hood with Air Flow
Tracking



In this Sequence:

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #4: Single Walk-in or Distillation Style Fume Hood with Air Flow Tracking

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has built-in local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the space pressure primary sensor method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an equal percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position.

If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for use by the laboratory pressure controller in calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Sequence of Operation



Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):
In addition, the controller provides the capability to


This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a

| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module.

remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.

15




NOTE:
SEE DRAWING A950054
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

WALK-IN OR DISTILLATION
TYPE FUME HOOD W/EXHAUST
FLOW MEASUREMENT

| | |
|---------------------------------------|--|
| DRAWING NO. A940044 | |
| SHEET 1 OF 1 | |
| DATE: 9-13-93 | |
| BY: [Signature] | |
| CHECKED: [Signature] | |
| APPROVED: [Signature] | |
| TYPICAL FUME HOOD LAB APPLICATION H 4 | |
| AUTO-FLOW | |
| SINGLE VALVE IN HOOD | |
| AIR FLOW TRACKING FLOW DIAGRAM | |

Application #5

Single Bench Fume Hood with
Auxiliary Makeup Supply Air



In this Sequence:

- *General*
- *Exhaust Control*
- *Auxiliary Makeup Air Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #5: Single Bench Fume Hood with Auxiliary Makeup Supply Air

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of displaying scaled point values and English language text point descriptors simultaneously
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Exhaust Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading acts as the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face

velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust

damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Auxiliary Makeup Air Control

The controller measures the exhaust air flow volume leaving the fume hood with a parallel plate pitot array in the fume hood exhaust duct. The controller also measures the auxiliary supply air being delivered to the hood with a parallel plate pitot array in the fume hood auxiliary supply air duct. These measurements are scaled into cubic feet per minute. The controller executes a process loop using a proportional, integral, derivative equation. The auxiliary supply air flow reading acts as the measured variable for the controller's PID loop. The setpoint for auxiliary supply air is seventy percent (adjustable) of the measured fume hood exhaust volume. The controller's PID output is characterized within the controller software in order to calculate an equal percentage control output to the fume hood auxiliary supply air damper. The fume hood auxiliary supply air damper is positioned as required by the controller with a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. If the supply air flow falls below the desired setpoint the controller modulates the normally closed fume hood auxiliary supply air damper toward the open position. If the supply air flow rises above the desired setpoint the controller modulates the normally closed fume hood auxiliary supply air damper toward the closed position. The fume hood exhaust flow measurement is also used by the laboratory pressure controller to calculate total lab exhaust volume (if required). The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has

Sequence of Operation



an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming


The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.

| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

The controller also allows for alarms on the exhaust and auxiliary supply air flow values based on user-adjustable limits. These alarms do not annunciate locally, they shall be transmitted over the communications network to the remote operator interface device.

[illegible]



NOTE:
SEE DRAWING PACKAGE A950055
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.


| | | | |
|---------------|--|-----------------------|--|
| DRAWN: LIP | | AUTO-FLOW | |
| ENGR: RVT | | REV: 10/01/01 | |
| CHKD: RVT | | SINGLE AIR BENCH HOOD | |
| DATE: 9-12-99 | | AIR FLOW DIAGRAM | |
| SHEET 1 | | DRAWING NO. A940045 | |

TYPICAL FUME HOOD LAB
APPLICATION # 5

BENCH TYPE FUME HOOD
W/AUXILIARY SUPPLY MAKE-UP AIR

Application #6

Single Walk-in or Distillation
Fume Hood with Auxiliary
Makeup Supply Air



In this Sequence:

- *General*
- *Exhaust Control*
- *Auxiliary Makeup Air Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #6: Single Walk-in or Distillation Fume Hood with Auxiliary Makeup Supply Air

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has built-in local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of displaying scaled point values and English language text point descriptors simultaneously
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Exhaust Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the space pressure primary sensor method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position.

If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Auxiliary Makeup Air Control

The controller measures the exhaust air flow volume leaving the fume hood with a parallel plate pitot array in the fume hood exhaust duct. The controller also measures the auxiliary supply air being delivered to the hood with a parallel plate pitot array in the fume hood auxiliary supply air duct. These measurements are scaled into cubic feet per minute. The controller executes a process loop using a proportional, integral, derivative equation. The auxiliary supply air flow reading acts as the measured variable for the controller's PID loop. The setpoint for auxiliary supply air is seventy percent (adjustable) of the measured fume hood exhaust volume. The controller's PID output is characterized within the controller software in order to calculate an equal percentage control output to the fume hood auxiliary supply air damper. The fume hood auxiliary supply air damper is positioned as required by the controller with a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. If the supply air flow falls below the desired setpoint the controller modulates the normally closed fume hood auxiliary supply air damper toward the open position. If the supply air flow rises above the desired setpoint the controller modulates the normally closed fume hood auxiliary supply air damper toward the closed position. The fume hood exhaust flow measurement is also used by the laboratory pressure controller to calculate total lab exhaust volume (if required). The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

Sequence of Operation



Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming


The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also allows for alarms on the exhaust and auxiliary supply air flow values based on user-adjustable limits. These alarms do not annunciate locally, they are transmitted over the communications network to the remote operator interface device.

| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

SINGLE WALK-IN HOOD w/Auxiliary Supply Air

[illegible]




NOTE:
SEE DRAWING A950056
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | |
|-----------------------------------------------------------------------|---------------|
| DRAWING NO. A940046 | |
| ENGINEER: RVT | DATE: 1-12-95 |
| CHECKED: RVT | SHEET: 1 |
| TYPICAL FUME HOOD LAB APPLICATION # 6 | |
| WALK-IN OR DISTILLATION TYPE FUME HOOD W/AUXILIARY SUPPLY MAKE-UP AIR | |

Application #7

Single Supply for Pressure
Control



In this Sequence:

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #7: Single Supply for Pressure Control

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels which supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network will have no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, and is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory using a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. If the measured supply air flow volume falls below its desired minimum air flow rate, as required to maintain the


necessary air changes per hour, the pressure control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume go above the desired maximum air flow rate the pressure control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Alarming

The controller provides the capability to alarm based on actuator failures related to either loss of power or actuator position feedback that does not meet commanded position. The controller also allows for alarms involving the supply air flow or total pressure differential measurements that are based on user adjustable limits. These alarms are not annunciated locally, they are transmitted via the communications network to the remote operator interface device.

26


NOTE:
SEE DRAWING A9S0059
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.



| | | | |
|-----------------------|-----------|-------------------------------|----|
| TYPICAL FUME HOOD LAD | | AUTID-ELDV | |
| APPLICATION # 7 | | SINGLE SUPPLY | |
| | | PRESSURE CONTROL FLOW DIAGRAM | |
| | | DRAWING NO. A940047 | |
| DESIGN LIP | ENGR. REV | DATE | BY |
| ENGR. REV | CHKD. REV | DATE | BY |
| ENGR. REV | CHKD. REV | DATE | BY |
| ENGR. REV | CHKD. REV | DATE | BY |

Application #8

Single Supply for Air Flow
Control



In this Sequence:

- *General*
- *Air Flow Control*
- *Volume Limit Control*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #8: Single Supply for Air Flow Control

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions independently of the network for all lab control aspects so that a failure of the communications network has no effect over laboratory pressure control.

Air Flow Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller executes a process loop using proportional, integral, and derivative control action. The supply air flow is the measured variable for the controller's PID loop. The setpoint for supply air flow is reset from local fume hood controllers or remotely from the BAS system. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned as required by the controller by a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop a minimum of five times per second. If the supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If the supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control


If the measured supply air flow volume falls below its desired minimum rate (adjustable) as required to maintain the necessary air changes per hour, the

supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. If the measured supply air flow volume rises above its desired maximum air flow rate (adjustable) the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Alarming

The controller provides the capability to alarm based on actuator failures related to either loss of power or actuator position feedback that does not meet commanded position. The controller also allows for alarms involving the supply air flow or total pressure differential measurements that are based on user adjustable limits. These alarms are not annunciated locally, they are transmitted via the communications network to the remote operator interface device.

[illegible]




NOTE:
SEE DRAWING A9500510
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|---------------------------------------|---------------|--------------------------------|--|
| DRAWING NO. A94004B | | AUTO-FLOW | |
| TYPICAL FUME HOOD LAB APPLICATION # 8 | | SINGLE SUPPLY | |
| DATE: 9-17-95 | | AIR FLOW TRACKING FLOW DIAGRAM | |
| SHEET | | DRAWING NO. A94004B | |
| DRAWN: LIP | ENGR: RWT | SHEET NO. 1 OF 1 | |
| CHECK: RWT | DATE: 9-17-95 | SHEET NO. 1 OF 1 | |

Application #9

Single Supply/General
Exhaust for Pressure Control
with Temperature Override



In this Sequence:

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Temperature Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW®

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #9: Single Supply/General Exhaust for Pressure Control with Temperature Override

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop a minimum of five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller executes a process loop using proportional, integral,

and derivative control action. The supply air flow reading is the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes a complete control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Temperature Override

The controller accepts an input from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired laboratory temperature. As the supply override signal increases, the controller increases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases, the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate).

Alarming


The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for alarms on the supply air flow and general exhaust air flow measurements based on

Sequence of Operation



user-adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

33




NOTE:
SEE DRAWINGS A9500513 & A9500515
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|-----------------------|-----------|-------------------------------|---------|
| TYPICAL FUME HOOD LAB | | AUTO-FLUX | |
| APPLICATION # 9 | | SINGLE SUPPLY/EXHAUST | |
| | | PRESSURE CONTROL FLOW DIAGRAM | |
| | | DRAWING NO. A940049 | |
| DRAWN: LIP | ENGR: RVT | CHKD: RVT | SHEET 1 |
| DATE: 9-13-95 | | | |

Application #10

Single Supply/General
Exhaust for Air Flow Control
with Temperature Override



In this Sequence:

- *General*
- *Air Flow Control*
- *Temperature Override*
- *Volume Limit Control*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #10: Single Supply/General Exhaust for Air Flow Control with Temperature Override

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Air Flow Control

The controller measures the laboratory supply air flow with a parallel plate pitot array in the laboratory supply air duct and it measures the exhaust air flow with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute (CFM). The controller executes its process loops using proportional, integral, and derivative control action. The supply air flow is the measured variable for one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal. The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value that is user adjustable. The exhaust air flow setpoint is based on the following equation:

$$\text{Exhaust CFM Setpoint} = \text{Current Supply CFM} + \Delta\text{CFM Setpoint}$$

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper

actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. When general exhaust air flow rises above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Temperature Override

The controller accepts an input from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired laboratory temperature. As the supply override signal increases, the controller increases the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases, the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its minimum value (the minimum design CFM for the required air change rate).

Volume Limit Control

Should the measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable), the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.


Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for alarms on the supply air flow and general exhaust air flow measurements based on user-adjustable limits. These alarms are not


Sequence of Operation



annunciated locally. They are transmitted over the communications network to the remote operator interface device.



NOTE:
SEE DRAWINGS A9500514 & A9500516
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.




| | | | |
|---------------|--|-------------------------------|--|
| DRAWN: LIP | | AUTO-FLOW | |
| ENGR: RVT | | SINGLE SUPPLY/EXHAUST | |
| CHECK: RVT | | AIR FLOW TRACING FLOW DIAGRAM | |
| DATE: 9-13-95 | | DRAWING NO. A9100410 | |
| SHEET 1 | | | |

TYPICAL FUME HOOD LAB
APPLICATION # 10

Application #11

Single Supply/General
Exhaust for Pressure Control
with Full Temperature Control



In this Sequence:

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Temperature Control*
- *Emergency Mode*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #11: Single Supply/General Exhaust for Pressure Control with Full Temperature Control

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network will have no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is -0.005"wc with respect to the adjacent corridor, which is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. When total pressure differential rises above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller executes a process loop using proportional, integral,

and derivative control action. The supply air flow reading is the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes a control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

$$\text{BTUs} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \times 1.085$$

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop controlling the general exhaust air damper. If the BTUs fall below the

Sequence of Operation



required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply).

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure, from any of the electronic actuators connected to the controller, related to either loss of power or actuator position feedback that does not meet the commanded position. The controller also allows for alarms on the supply air flow, general exhaust air flow, supply duct temperature, space temperature or total pressure differential measurements based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted over the communications network to the remote operator interface device.




NOTE:
SEE DRAWINGS A9500511 & A9500515
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|-----------------------|--|----------------------------|--|
| TYPICAL FUME HOOD LAB | | AUTO-FLOW | |
| APPLICATION # 11 | | BY J. H. HARRIS | |
| | | ENGINEER, RVT | |
| | | SINGLE SUPPLY/EXHAUST | |
| | | PRESSURE/TEMP FLOW DIAGRAM | |
| | | DRAWING NO. A9400411 | |

Application #12

Single Supply/General
Exhaust for Air Flow Control
with Full Temperature Control



In this Sequence:

- *General*
- *Air Flow Control*
- *Temperature Control*
- *Volume Limit Control*
- *Emergency Mode*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW®

Sequence of Operation

Application #12: Single Supply/General Exhaust for Air Flow Control with Full Temperature Control

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels which supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Air Flow Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct and measures the exhaust air flow leaving the space with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute. The controller executes its process loops with proportional, integral, and derivative control action. The supply air flow is the measured variable for one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature control. The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value which is user adjustable. The exhaust air flow setpoint is based on the following equation:

$$\text{Exhaust CFM Setpoint} = (\text{Current Supply CFM} + \Delta \text{CFM Setpoint})$$

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller

Laboratory and Fume Hood Controls Engineering Guide

through an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. If general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. On a rise in general exhaust air flow above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures, combined with the measured supply air flow rate, are used to calculate the instantaneous BTUs used by the space to maintain temperature. The equation for this process is as follows

$$\text{BTUs} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \times 1.085$$

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop to control the laboratory supply air damper. If BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop, controlling the supply air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint, the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop, controlling the supply air damper toward its

Sequence of Operation



maximum value (the maximum design CFM for the lab supply).

Volume Limit Control

Should the measured supply air flow volume fall below its desired minimum air flow rate as required to maintain the necessary air changes per hour, the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. If the measured supply air flow volume rises above its desired maximum air flow rate the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Emergency Mode


The controller has the capability to accept an external contact from an emergency push button to indicate the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure from any of the electronic actuators connected to it that is related to either loss of power or to actuator position feedback that does not meet commanded position. Additionally, the controller also allows for alarms related to the supply air flow, general exhaust air flow, supply duct temperature and space temperature measurements. The setpoints for these are alarms are based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted over the communications network to the remote operator interface device.

Application #13

Single Supply for Pressure
Control with Single Bench
Fume Hood



In this Sequence:

LABORATORY

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Alarming*

FUME HOOD

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW®

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #13: Single Supply for Pressure Control with Single Bench Fume Hood

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading acts as the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory using a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. If the measured supply air flow volume falls below its desired minimum air flow rate, as required to maintain the necessary air changes per hour, the pressure control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume go above the desired maximum air flow rate the pressure control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Alarming

The controller provides the capability to alarm based on actuator failures related to either loss of power or actuator position feedback that does not meet commanded position. The controller also allows for alarms involving the supply air flow or total pressure differential measurements that are based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted via the communications network to the remote operator interface device.

FUME HOOD

General

The fume hood is equipped with a stand alone networkable microprocessor- based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key

Sequence of Operation

- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, which is adjustable either locally or

| | | | |
|----------------|------------------------|----------------------------------|---------------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

remotely. The controller's PID output is calculated by the controller software in order to calculate an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller, after an adjustable time delay, generates a local alarm as described in the alarming section.


Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1) :

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.




NOTE:
SEE DRAWINGS A950051 & A950059
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|----------------------------------------|--|----------------------|----------------|
| TYPICAL FUME HOOD LAB APPLICATION # 13 | | DESIGN LIP | AUTO-FLOW |
| | | ENGR. BY | ENGR. BY |
| | | CHECK BY | CHECK BY |
| | | DATE: 12-15-15 | DATE: 12-15-15 |
| | | SHEET | SHEET |
| | | DRAWING NO. A9400413 | |

BENCH TYPE FUME HOOD

Application #14

Single Supply for Air Flow
Control with Single Bench
Fume Hood



In this Sequence:

LABORATORY

- *General*
- *Air Flow Control*
- *Volume Limit Control*
- *Alarming*

FUME HOOD

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #14: Single Supply for Air Flow Control with Single Bench Fume Hood

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions independently of the network for all lab control aspects so that a failure of the communications network has no effect over laboratory pressure control.

Air Flow Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into units of cubic feet per minute. The controller executes a process loop using proportional, integral, and derivative control action. The supply air flow is the measured variable for the controller's PID loop. The setpoint for supply air flow is reset to satisfy the fume hood exhaust volume as read from the cascaded output of the fume hood controller minus a set user-adjustable differential volume. The controller's PID output is used to control the laboratory supply damper. The laboratory supply damper is positioned as required by the controller by an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop of five times per second. If the supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If the supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

If the measured supply air flow volume falls below its desired minimum rate (adjustable) as required to maintain the necessary air changes per hour, the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. If the measured supply air flow volume rises above its desired maximum air flow rate (adjustable) the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Alarming

The controller provides the capability to alarm based on actuator failures related to either loss of power or actuator position feedback that does not meet commanded position. The controller also allows for alarms involving the supply air flow or total pressure differential measurements that are based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted via the communications network to the remote operator interface device.

FUME HOOD

General

The fume hood is equipped with a stand alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all

Sequence of Operation



hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set at 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its

totalized signal to either another fume hood controller or to the laboratory pressure controller.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | |
|------------------------|---------------------|--------------------------------------|
| Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |


Table 1

control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood via a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute (CFM) for use by the laboratory pressure controller in calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value that is based on user-adjustable limits. This alarm is not annunciated locally, but is transmitted over the communications network to the remote operator interface device.




NOTE:
SEE DRAWINGS A950053 & A9500510
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

BENCH TYPE FUME HOOD

| | | | |
|-----------------------|-----------|---------------------------------|----------|
| TYPICAL FUME HOOD LAB | | AUTO-FLOW | |
| DESIGN: LVP | CHKD: RVT | DATE: 8-13-95 | SHEET: 1 |
| APPLICATION H 14 | | DRAWING NO. A94000414 | |
| | | SINGLE SUPPLY/SINGLE BENCH HOOD | |
| | | AIR FLOW TRACKING FLOW DIAGRAM | |

Application #15

Single Supply/General
Exhaust for Pressure Control
with Temperature Override
and Single Bench Fume Hood



In this Sequence:

LABORATORY

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Temperature Override*
- *Emergency Mode*
- *Alarming*

FUME HOOD

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #15: Single Supply/General Exhaust for Pressure Control with Temperature Override and Single Bench Fume Hood

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop of five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller executes a process loop using proportional, integral, and derivative control action. The supply air flow reading is the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Temperature Override

The controller accepts an input from the Building Automation System (BAS) to determine the amount of supply flow override required to maintain the desired temperature in the laboratory. As the supply override signal increases the controller increases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate).

Sequence of Operation



Emergency Mode

The controller has the capability to accept an external contact from an emergency push button used to indicate the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure, from any of the electronic actuators connected to the controller, related to either loss of power or actuator position feedback that does not meet the commanded position. The controller also allows for alarms on the supply air flow, general exhaust air flow, supply duct temperature, space temperature or total pressure differential measurements based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted over the communications network to the remote operator interface device.

FUME HOOD

General

The fume hood is equipped with a stand alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the

communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop of twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming


The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a

| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by

remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.




BENCH TYPE FUME HOOD

| | |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| TYPICAL FUME HOOD LAB APPLICATION # 15 | DRAWING NO. A94000415 SHEET 15 DATE: 9-11-93 CHECKED BY DESIGNED BY ENGINEER BY AUTO-CAD FILED BY DRAWN BY DATE: 9-11-93 SHEET 15 DRAWING NO. A94000415 |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

Application #16

Single Supply/General
Exhaust for Air Flow Control
with Temperature Override
and Single Bench Hood



In this Sequence:

LABORATORY

- *General*
- *Air Flow Control*
- *Temperature Override*
- *Volume Limit Control*
- *Emergency Mode*
- *Alarming*

FUME HOOD

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #16: Single Supply/General Exhaust for Air Flow Control with Temperature Override and Single Bench Hood

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Air Flow Control

The controller measures the laboratory supply air flow with a parallel plate pitot array in the laboratory supply air duct and it measures the exhaust air flow with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute (CFM). The controller executes its process loops using proportional, integral, and derivative control action. The supply air flow acts as the measured variable for the one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal or the air flow rate required to satisfy the fume hood makeup air volume, whichever is greater. The supply makeup air flow rate is determined by the following equation:

Supply Make-Up Flow Setpoint CFM = (Total Hood Exhaust Air Flow CFM + Minimum General Exhaust Air Flow CFM - Desired Differential CFM Setpoint)

The controller's PID output is used to control output to the laboratory supply air damper. The laboratory supply air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for

positive positioning. The controller executes its control loop of five times per second. When supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the open position. When supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the closed position.

The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value that is user adjustable. The exhaust air flow setpoint is based on the following equation:

Exhaust CFM Setpoint = Current Supply CFM + Δ CFM Setpoint - Total Hood Exhaust CFM

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop of five times per second. When general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. When general exhaust air flow rises above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Temperature Override

The controller accepts an input from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired laboratory temperature. As the supply override signal increases, the controller increases the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases, the controller decreases the minimum supply air flow setpoint associated with

Sequence of Operation

the PID loop controlling the supply air damper toward its minimum value (the minimum design CFM for the required air change rate).

Volume Limit Control

Should the measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable), the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for the alarming of the supply air flow and general exhaust air flow measurements based on user-adjustable limits. These alarms are not annunciated locally. They are transmitted via the communications network to the remote operator interface device.

FUME HOOD

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local



operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of visually displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop of twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1) :


| | | | |
|----------------|------------------------|----------------------------------|---------------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

use by the laboratory pressure controller in calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.




BENCH TYPE FUME HOOD

| | | | |
|-----------------------|---------------|---------------------------------|--|
| TYPICAL FUME HOOD LAB | | AUTO-CLAVE | |
| APPLICATION # 16 | | SINGLE SUPPLY/EXHAUST FUME HOOD | |
| SHAW-WALKER | ENGRD BY | AIR FLOW | |
| CHGR BY | DATE: 9-12-95 | FLOW RANGE | |
| SHEET | DRAWING NO. | A9400416 | |

Application #17

Single Supply/General
Exhaust for Pressure Control
with Full Temperature Control
and Single Bench Fume Hood



In this Sequence:

LABORATORY

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Temperature Control*
- *Emergency Mode*
- *Alarming*

FUME HOOD

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #17: Single Supply/General Exhaust for Pressure Control with Full Temperature Control and Single Bench Fume Hood

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller

executes a process loop using proportional, integral, and derivative control action. The supply air flow reading acts as the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

$$\text{BTU's} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \times 1.085$$

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop controlling the general exhaust air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow

setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate).

As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply).

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on actuator failures related to either loss of power or actuator position feedback that does not meet commanded position. The controller also allows for alarms involving the supply air flow or total pressure differential measurements that are based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted via the communications network to the remote operator interface device.

FUME HOOD

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors

Sequence of Operation

- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In this mode the controller displays audible alarm sounds and the

Sequence of Operation





red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.




NOTE:
SEE DRAWINGS A950051, A9500511 &
A9500515 FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | | | | | |
|-----------------------|--|------------------------------|--|----------------------|--|----------------------|--|
| DRAWN LIP | | CHECKED RVT | | DATE: 9-11-95 | | SHEET | |
| ENGINEER RVT | | CHECKED RVT | | DATE: 9-11-95 | | SHEET | |
| TYPICAL FUME HOOD LAB | | APPLICATION # 17 | | DRAWING NO. A9400417 | | AUTO-FLD | |
| BY: 10/01/95 | | SHEET SUPPLY/EXHAUST CONTROL | | PRESSURE/TEMPERATURE | | DRAWING NO. A9400417 | |

BENCH TYPE FUME HOOD

Application #18

Single Supply/General
Exhaust for Air Flow Control
with Full Temperature Control
and Single Bench Hood



In this Sequence:

LABORATORY

- *General*
- *Air Flow Control*
- *Temperature Control*
- *Volume Limit Control*
- *Emergency Mode*
- *Alarming*

FUME HOOD

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #18: Single Supply/General Exhaust for Air Flow Control with Full Temperature Control and Single Bench Hood

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Air Flow Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct and measures the exhaust air flow leaving the space with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute. The controller executes its process loops with proportional, integral, and derivative control action. The supply air flow is the measured variable for the one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal or the air flow rate required to satisfy the fume hood makeup air volume, whichever is greater. The supply makeup air flow rate is determined by the following equation:

Supply Makeup Flow Setpoint CFM = (Total Hood Exhaust Air Flow CFM + Minimum General Exhaust Air Flow CFM - Desired Differential CFM Setpoint)

The controller's PID output is used to control output to the laboratory supply air damper. The laboratory supply air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When supply air

flow falls below the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the open position. When supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the closed position.

The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value which is adjustable by the user. The exhaust air flow setpoint is based on the following equation:

Exhaust CFM Setpoint = (Current Supply CFM + Differential CFM Setpoint - Total Hood Exhaust CFM)

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller through an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. If general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. On a rise in general exhaust air flow above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

BTU's = [(Supply Discharge Temperature - Space Temperature) × Supply Air Flow Rate] × 1.085

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired

Sequence of Operation

instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop controlling the general exhaust air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the laboratory supply air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its maximum value (the maximum design CFM for the lab supply).

Volume Limit Control

Should the measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable), the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming



The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for the alarming of the supply air flow and general exhaust air flow measurements based on user adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

FUME HOOD

General

The fume hood is equipped with a stand alone networkable microprocessor- based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The

damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm

Table 1 (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for use by the laboratory pressure controller in calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.


Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the red controller display's LED and the

audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.

| | | |
|------------------------|---------------------|--------------------------------------|
| Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |



NOTE:
SEE DRAWINGS A950053, A9500512 &
A9500516 FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.


BENCH TYPE FUME HOOD

| | |
|-------------|------------------|
| DESIGN LIP | AUTO-FLOW |
| ENGINEER | W. J. YOUNG, JR. |
| CHECKED | W. J. YOUNG, JR. |
| DATE | 5-13-75 |
| SHEET | 1 |
| PROJECT | FLUOR DUBOIS |
| DRAWING NO. | A940041B |

TYPICAL FUME HOOD LAB
APPLICATION # 18

Application #19

Single Supply/General
Exhaust for Pressure Control
with Temperature Override
and Dual Bench Fume Hoods



In this Sequence:

LABORATORY

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Temperature Override*
- *Emergency Mode*
- *Alarming*

FUME HOODS

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #19: Single Supply/General Exhaust for Pressure Control with Temperature Override and Dual Bench Fume Hoods

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect over laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller

executes a process loop using proportional, integral, derivative control action. The supply air flow reading is the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes a control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Temperature Override

The controller accepts an input from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired laboratory temperature. As the supply override signal increases, the controller increases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases, the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate).

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications

network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure, from any of the electronic actuators connected to the controller, related to either loss of power or actuator position feedback that does not meet the commanded position. The controller also allows for alarms on the supply air flow or total pressure differential measurements based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted over the communications network to the remote operator interface device.

FUME HOODS

General

The fume hood is equipped with a stand alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Sequence of Operation

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, which is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an equal percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Sequence of Operation



Alarming


The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1) :

audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller’s display module. This key does not cancel the alarm condition but does silence the alarm for a time period. After the time delay the audible alarm returns if the condition

| | | | |
|---------|-----------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display’s red LED and the

has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.




NOTE:
SEE DRAWINGS A950051, A9500513 &
A9500515 FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|-----------------------|------------|------------------|---------|
| TYPICAL FUME HOOD LAB | | APPLICATION # 19 | |
| DRAWN LIP | ENG'D. REV | DATE 9-11-90 | SHEET 1 |
| CHECKED BY | | DATE 9-11-90 | |
| APPROVED BY | | DATE 9-11-90 | |
| DRAWING NO. | | A9400419 | |

2 BENCH TYPE FUME HOODS

Application #20

Single Supply/General
Exhaust for Air Flow Control
with Temperature Override
and Dual Bench Hoods



In this Sequence:

LABORATORY

- *General*
- *Air Flow Control*
- *Temperature Override*
- *Volume Limit Control*
- *Emergency Mode*
- *Alarming*

FUME HOODS

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #20: Single Supply/General Exhaust for Air Flow Control with Temperature Override and Dual Bench Hoods

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Air Flow Control

The controller measures the laboratory supply air flow with a parallel plate pitot array in the laboratory supply air duct and it measures the exhaust air flow with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute (CFM). The controller executes its process loops using proportional, integral, and derivative control action. The supply air flow acts as the measured variable for the one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal or the air flow rate required to satisfy the fume hood makeup air volume, whichever is greater. The supply makeup air flow rate is determined by the following equation:

Supply Make-Up Flow Setpoint CFM = (Total Hood Exhaust Air Flow CFM + Minimum General Exhaust Air Flow CFM - Desired Differential CFM Setpoint)

The controller's PID output is used to control output to the laboratory supply air damper. The laboratory supply air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for

positive positioning. The controller executes its control loop five times per second. When supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the open position. When supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the closed position.

The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value that is user adjustable. The exhaust air flow setpoint is based on the following equation:

Exhaust CFM Setpoint = Current Supply CFM + Δ CFM Setpoint - Total Hood Exhaust CFM

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop five times per second. When general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. When general exhaust air flow rises above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Temperature Override

The controller accepts an input from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired laboratory temperature. As the supply override signal increases, the controller increases the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases, the controller decreases the

Sequence of Operation

minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its minimum value (the minimum design CFM for the required air change rate).

Volume Limit Control

Should the measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable), the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for the alarming of the supply air flow and general exhaust air flow measurements based on user adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

FUME HOODS

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local

*Laboratory and Fume Hood
Controls Engineering Guide*

operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for use by the laboratory pressure controller in

Alarming

The controller has the ability to initiate alarms both locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|----------------|------------------------|----------------------------------|---------------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |


calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user-adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.

84



NOTE:
SEE DRAWINGS A950053, A9500514 &
A9500516 FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.


2 BENCH TYPE FUME HOODS

| | | |
|----------------------|-------------|---------------------------|
| DESIGN LIP | AUTO-FLOW | EX 10000000 & EX 10000000 |
| ENGR: RVT | EX 10000000 | EX 10000000 |
| CHECK: RVT | EX 10000000 | EX 10000000 |
| DATE: 9-13-94 | EX 10000000 | EX 10000000 |
| SHEET | EX 10000000 | EX 10000000 |
| DRAWING NO. A9400420 | | |

TYPICAL FUME HOOD LAB
APPLICATION # 20

Application #21

Single Supply/General
Exhaust for Pressure Control
with Full Temperature Control
and Dual Bench Fume Hoods



In this Sequence:

LABORATORY

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Temperature Control*
- *Emergency Mode*
- *Alarming*

FUME HOODS

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #21: Single Supply/General Exhaust for Pressure Control with Full Temperature Control and Dual Bench Fume Hoods

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop of five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller executes a process loop using proportional, integral,

and derivative control action. The supply air flow reading is the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes a complete control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

$$\text{BTUs} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \times 1.085$$

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop to control the general exhaust air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the

general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply).

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for the alarming of the supply air flow and general exhaust air flow measurements based on user-adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

FUME HOODS

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors

Sequence of Operation

- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, which is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an equal percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Sequence of Operation



Emergency Override


The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1) :

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.

| | | | |
|---------|-----------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |




NOTE:
SEE DRAWINGS A950051, A9500511 &
A9500515 FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|-----------------------|----------|----------------------|--|
| TYPICAL FUME HOOD LAB | | DRAWING NO. A9400421 | |
| APPLICATION # 21 | | SHEET 1 | |
| DESIGN: LIP | ENG: RY | DATE: 9-13-95 | |
| CHKD: RY | CHKD: RY | | |
| BY: J. L. LIPMAN | | BY: J. L. LIPMAN | |
| DATE: 9-13-95 | | DATE: 9-13-95 | |
| SHEET 1 | | SHEET 1 | |

2 BENCH TYPE FUME HOODS

Application #22

Single Supply/General
Exhaust for Air Flow Control
with Full Temperature Control
and Dual Bench Hoods



In this Sequence:

LABORATORY

- *General*
- *Air Flow Control*
- *Temperature Control*
- *Volume Limit Control*
- *Emergency Mode*
- *Alarming*

FUME HOODS

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW®

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #22: Single Supply/General Exhaust for Air Flow Control with Full Temperature Control and Dual Bench Hoods

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Air Flow Control

The controller measures the laboratory supply air flow with a parallel plate pitot array in the laboratory supply air duct and it measures the exhaust air flow with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute (CFM). The controller executes its process loops using proportional, integral, and derivative control action. The supply air flow is the measured variable for the one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal or the air flow rate required to satisfy the fume hood makeup air volume, whichever is greater. The supply makeup air flow rate is determined by the following equation:

Supply Make-Up Flow Setpoint CFM = (Total Hood Exhaust Air Flow CFM + Minimum General Exhaust Air Flow CFM - Desired Differential CFM Setpoint)

The controller's PID output is used to control output to the laboratory supply air damper. The laboratory supply air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for

positive positioning. The controller executes its control loop of five times per second. When supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the open position. When supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory supply air damper toward the closed position.

The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value that is user adjustable. The exhaust air flow setpoint is based on the following equation:

Exhaust CFM Setpoint = Current Supply CFM + Δ CFM Setpoint-Total Hood Exhaust CFM

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. When general exhaust air flow rises above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

BTUs = [(Supply Discharge Temperature - Space Temperature) \times Supply Air Flow Rate] \div 1.085

Sequence of Operation

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop controlling the laboratory supply air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its maximum value (the maximum design CFM for the lab supply).

Volume Limit Control

Should the measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable), the supply flow control PID loop is overridden and the laboratory supply damper is controlled to maintain the desired maximum air flow.

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this

position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for the alarming of the supply air flow and general exhaust air flow measurements based on user-adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

FUME HOODS

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect over fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the

measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, which is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms both locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|----------------|------------------------|----------------------------------|---------------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |


normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for use by the laboratory pressure controller in calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user-adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.

PROJECT NO.:




NOTE:
SEE DRAWINGS A950053, A9500512 &
A9500516 FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|-----------------------|----------|----------------------|-------|
| TYPICAL FUME HOOD LAB | | DRAWING NO. A9400422 | |
| DESIGN | ENGINEER | DATE | SHEET |
| CHECKED | DATE | DATE | SHEET |
| APPLICATION # 22 | | DRAWING NO. A9400422 | |

2 BENCH TYPE FUME HOODS

Application #23

Single Supply/General
Exhaust for Pressure Control
with Temperature Override
and Pressure Controlled
Internal Office



In this Sequence:

- *General*
- *Lab Pressure Control*
- *Lab Office Pressure Control*
- *Lab Volume Limit Control*
- *Lab Office Volume Limit Control*
- *Lab Temperature Override*
- *Lab Office Temperature Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #23: Single Supply/General Exhaust for Pressure Control with Temperature Override and Pressure Controlled Internal Lab Office

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Lab Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading acts as the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Lab Office Pressure Control

The total pressure differential is measured between the laboratory office space and the surrounding laboratory. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the

measured variable for the controller's PID loop. The setpoint for total pressure differential is + 0.005"wc with respect to the surrounding laboratory, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory office supply damper. The laboratory office supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory office supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory office supply damper toward the closed position.

Lab Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller executes a process loop using proportional, integral, and derivative control action. The supply air flow reading is the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Lab Office Volume Limit Control

The controller measures the supply air flow being delivered to the laboratory office using a parallel plate pitot array in the laboratory office supply air duct. This measurement is scaled into cubic feet per minute. The setpoint for supply air flow rate is that which maintains the desired minimum air changes per hour. The setpoint is adjustable either locally or remotely. The controller resets the total pressure differential PID loop setpoint in order to maintain the minimum or maximum supply flow limits. When the supply air flow rate falls below the desired minimum setpoint the controller resets the desired differential setpoint upward in order to increase supply flow. When the supply air flow rate rises above the desired minimum setpoint the controller resets the desired setpoint downward in order to decrease supply air flow.

Lab Temperature Override

The controller accepts an input from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired laboratory temperature. As the supply override signal increases, the controller increases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases, the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate).

Lab Office Temperature Override


The controller accepts an input from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired laboratory office temperature. As the supply override signal increases the controller increases the minimum supply air flow setpoint associated with the reset of the differential setpoint. As the supply override signal decreases the controller decreases the minimum supply air flow setpoint associated with the reset of the differential setpoint.

Alarming

The controller provides the capability to alarm based on an actuator failure, related to either loss of power or actuator position feedback not meeting its commanded position. Additionally the controller also allows for the alarming of either of the supply air flows or total pressure differential measurements based on user-adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

Application #24

Single Supply/General
Exhaust for Air Flow Control
with Temperature Override
and Flow Controlled Internal
Lab Office



In this Sequence:

- *General*
- *Lab Air Flow Control*
- *Lab Office Air Flow Control*
- *Temperature Override*
- *Volume Limit Control*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #24: Single Supply/General Exhaust for Air Flow Control with Temperature Override and Flow Controlled Internal Lab Office

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Lab Air Flow Control

The controller measures the laboratory supply air flow with a parallel plate pitot array in the laboratory supply air duct and it measures the exhaust air flow with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute (CFM). The controller executes its process loops using proportional, integral, and derivative control action. The supply air flow is the measured variable for one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal. The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value that is user-adjustable. The exhaust air flow setpoint is based on the following equation:

$$\text{Exhaust CFM Setpoint} = \text{Current Supply CFM} + \text{Current Office Supply CFM} + \Delta\text{CFM Setpoint}$$

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller

with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. When general exhaust air flow rises above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Lab Office Air Flow Control

The controller measures the supply air flow delivered to the laboratory office with a parallel plate pitot array in the laboratory office supply air duct. This measurement is scaled into units of cubic feet per minute. The controller executes a process loop using proportional, integral, and derivative control action. The lab office supply air flow is the measured variable for the one of the controller's PID loops. The setpoint for lab office supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal. The damper actuator incorporates internal feedback for positive positioning. If the laboratory office supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory office supply air damper toward the open position. If the laboratory office supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory office supply air damper toward the closed position.

Temperature Override

The controller accepts inputs from the Building Automation System (BAS) which determines the amount of supply flow override required to maintain the desired temperature in the laboratory and the laboratory office. As either supply override signal increases the controller increases the respective minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its maximum value (the maximum design CFM for the lab supply). As the supply override signal decreases, the controller decreases the respective minimum supply air flow setpoint associated with the PID loop controlling the supply

air damper toward its minimum value (the minimum design CFM for the required air change rate).

Volume Limit Control


Should either laboratory or laboratory office measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, the respective supply flow control PID loop is overridden and the respective supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable) the respective supply flow control PID loop is overridden and the respective supply damper is controlled to maintain the desired maximum air flow.

Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for alarms on either of the supply air flows and the general exhaust air flow measurements based on user-adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

Application #25

Single Supply/General
Exhaust for Pressure Control
with Full Temperature Control
and Internal Lab Office with
Pressure Control and Full
Temperature Control



In this Sequence:

- *General*
- *Lab Pressure Control*
- *Lab Office Presssure Control*
- *Lab Volume Limit Control*
- *Lab Office Volume Limit Control*
- *Lab Temperature Control*
- *Lab Office Temperature Control*
- *Emergency Mode*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #25: Single Supply/General Exhaust for Pressure Control with Full Temperature Control and Internal Lab Office with Pressure Control and Full Temperature Control

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Lab Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Lab Office Pressure Control

The total pressure differential is measured between the laboratory office space and the surrounding laboratory. The controller executes a process loop using proportional, integral, and derivative control

action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is + 0.005"wc with respect to the surrounding laboratory, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory office supply damper. The laboratory office supply damper is positioned by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory office supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory office supply damper toward the closed position.

Lab Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in the laboratory supply air duct. This measurement is scaled into cubic feet per minute. The controller executes a process loop using proportional, integral, and derivative control action. The supply air flow reading is the measured variable for the controller's PID loop. The setpoint for supply air flow rate is that which maintains the desired minimum number of air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes a control loop five times per second. When the supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Lab Office Volume Limit Control

The controller measures the supply air flow being delivered to the laboratory office using a parallel plate pitot array in the laboratory office supply air duct. This measurement is scaled into cubic feet per minute. The setpoint for supply air flow rate is that which maintains the desired minimum air changes per hour. The setpoint is adjustable either locally or remotely. The controller resets the total pressure differential PID loop setpoint in order to maintain the minimum or maximum supply flow limits. When the supply air flow rate falls below the desired minimum setpoint the controller resets the desired differential setpoint upward in order to increase supply flow. When the supply air flow rate rises above the desired minimum setpoint the controller resets the desired setpoint downward in order to decrease supply air flow.

Lab Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

$$\text{BTUs} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \div 1.085$$

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop controlling the general exhaust air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the

controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply).

Lab Office Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory office space and the temperature of the laboratory office space itself. These temperatures combined with the measured laboratory office supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain temperature. The equation for this process is as follows:

$$\text{BTU's} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \div 1.085$$

One of the controllers PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum laboratory office supply flow setpoint and the associated PID loop. If the BTUs fall below the required setpoint the controller decreases the minimum laboratory office supply air flow toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum laboratory office supply air flow setpoint toward its maximum value (the maximum design CFM for the lab office supply).

Emergency Mode


The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-

Sequence of Operation

adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure, from any of the electronic actuators connected to the controller, related to either loss of power or actuator position feedback that does not meet the commanded position. The controller also allows for alarms on the supply air flow, general exhaust air flow, supply duct temperature, space temperature or total pressure differential measurements based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted over the communications network to the remote operator interface device.




NOTE:
SEE DRAWINGS A950051D & A950051S
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| TYPICAL FUME HOOD LAB | | | | AUTOMATIC FLOW | |
|-----------------------|---------|-----|------|----------------------|---------|
| REVISION | DATE | BY | CHKD | DATE | BY |
| 1 | 9-13-93 | SPJ | | 1 | 9-13-93 |
| APPLICATION # 25 | | | | DRAWING NO. A9400425 | |

Application #26

Single Supply/General
Exhaust for Air Flow Control
with Full Temperature Control
and Internal Lab Office with
Air Flow Control and Full
Temperature Control



In this Sequence:

- *General*
- *Lab Air Flow Control*
- *Lab Office Air Flow Control*
- *Lab Temperature Control*
- *Lab Office Temperature Control*
- *Volume Limit Control*
- *Emergency Mode*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #26: Single Supply/General Exhaust for Air Flow Control with Full Temperature Control and Internal Lab Office with Air Flow Control and Full Temperature Control

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network will have no effect on laboratory pressure control.

Lab Air Flow Control

The controller measures the laboratory supply air flow with a parallel plate pitot array in the laboratory supply air duct and it measures the exhaust air flow with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute (CFM). The controller executes its process loops using proportional, integral, and derivative control action. The supply air flow acts as the measured variable for one of the controller's PID loops. The setpoint for supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal. The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value that is user adjustable. The exhaust air flow setpoint is based on the following equation:

$$\text{Exhaust CFM Setpoint} = \text{Current Lab Supply CFM} + \text{Current Lab Office CFM} + \Delta\text{CFM Setpoint}$$

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller

with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop five times per second. When general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. When general exhaust air flow rises above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Lab Office Air Flow Control

The controller measures the supply air flow delivered to the laboratory office with a parallel plate pitot array in the laboratory office supply air duct. This measurement is scaled into units of cubic feet per minute. The controller executes a process loop using proportional, integral, and derivative control action. The lab office supply air flow is the measured variable for one of the controller's PID loops. The setpoint for lab office supply air flow is the minimum air flow rate required for the desired air changes per hour, reset from the temperature override signal. The damper actuator incorporates internal feedback for positive positioning. If the laboratory office supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory office supply air damper toward the open position. If the laboratory office supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory office supply air damper toward the closed position.

Lab Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

$$\text{BTUs} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \times 1.085$$

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is

Sequence of Operation

used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop controlling the laboratory supply air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the supply air damper toward its maximum value (the maximum design CFM for the lab supply).

Lab Office Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory office space and the temperature of the laboratory office space itself. These temperatures combined with the measured laboratory office supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

$$\text{BTUs} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \times 1.085.$$

One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coil. Additionally the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum laboratory office supply flow setpoint and the associated PID loop. If the BTUs fall below the required setpoint the controller decreases the minimum laboratory office supply air flow toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If the calculated

BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum laboratory office supply air flow setpoint toward its maximum value (the maximum design CFM for the lab office supply).

Volume Limit Control


Should the measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, the respective supply flow control PID loop is overridden and the respective supply damper is controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable), the respective supply flow control PID loop is overridden and the respective supply damper is controlled to maintain the desired maximum air flow.

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming


The controller provides the capability to alarm based on an actuator failure, from any of the electronic actuators connected to the controller, related to either loss of power or actuator position feedback that does not meet the commanded position. The controller also allows for alarms on the supply air flow, general exhaust air flow, supply duct temperature, space temperature or total pressure differential measurements based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted over the communications network to the remote operator interface device.



| | |
|----------------------------------------|--|
| DRAWING NO. A9400426 | |
| SHEET 1 | |
| DATE: 9-13-95 | |
| CHECKED BY: [Signature] | |
| DESIGNED BY: [Signature] | |
| TYPICAL FUME HOOD LAB APPLICATION # 26 | |
| AUTO-FLY | |
| AIR FLOW TRANSDUCING CONTROL | |
| BY: [Signature] | |
| DATE: 9-13-95 | |
| DRAWING NO. A9400426 | |

Application #27

Dual Supply Single General
Exhaust for Pressure Control
with Full Temperature Control
and Four Bench Fume Hoods



In this Sequence:

LABORATORY

- *General*
- *Pressure Control*
- *Volume Limit Control*
- *Temperature Control*
- *Emergency Mode*
- *Alarming*

FUME HOODS

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #27: Dual Supply Single General Exhaust for Pressure Control with Full Temperature Control and Four Bench Fume Hoods

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading acts as the measured variable for the controller's PID loop. The setpoint for total pressure differential is -0.005"wc with respect to the adjacent corridor, which is adjustable either locally or remotely. The controller's PID output is used to control each output to the laboratory supply dampers. The laboratory supply dampers are positioned as required by the controller by electronic damper actuators. The damper actuators incorporate internal feedback for positive positioning. The controller executes its control loop five times per second. When total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. When total pressure differential rises above the desired setpoint the controller modulates the normally closed laboratory supply dampers toward the closed position.

Volume Limit Control

The controller measures the supply air flow delivered to the laboratory with a parallel plate pitot array in each of the laboratory supply air ducts. This measurement is scaled into cubic feet per minute. The controller executes a process loop using a proportional, integral, and derivative control action. The total supply air flow reading acts as the measured variable for the controller's PID loop. The setpoint for total supply air flow rate shall be that which maintains the desired minimum air changes per hour. The setpoint is adjustable either locally or remotely. The controller's PID output is used to control the output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes a control loop five times per second. When the total supply air flow rate falls below the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the open position. Opening the damper causes a drop in total pressure differential that results in an increase in supply air delivered to the laboratory. When the total supply air flow rate goes above the desired minimum setpoint the controller modulates the normally open laboratory general exhaust damper toward the closed position. This action results in an increase in total pressure differential causing a reduction in the supply air delivered to the lab.

Temperature Control

The controller measures the temperature of the supply air in each of the supply ducts serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain temperature. The equation for this process is as follows:

$$\text{BTUs} = [(\text{Supply Discharge Temperature} - \text{Space Temperature}) \times \text{Supply Air Flow Rate}] \times 1.085$$

This equation is calculated for each supply duct serving the lab space and summed to provide a total BTU result. One of the controller's PID loops is used

in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. This BTU setpoint is used in conjunction with the calculated current BTUs and a PID loop to control the laboratory supply duct reheat coils. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and the associated PID loop controlling the general exhaust air damper. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the output to the reheat coil valve toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil output toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the PID loop controlling the general exhaust air damper toward its maximum value (the maximum design CFM for the lab supply).

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure, from any of the electronic actuators connected to the controller, related to either loss of power or actuator position feedback that does not meet the commanded position. The controller also allows for alarms on the supply air flow, supply duct temperature, space temperature or total pressure differential measurements based on user-adjustable limits. These alarms are not annunciated locally, they are transmitted over the communications network to the remote operator interface device.

FUME HOODS

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an

Sequence of Operation

adjustable time delay the controller generates a local alarm (as described in the alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|----------------|------------------------|----------------------------------|---------------------------------------------|
| Table 1 | Velocity Unsafe | <i><80 FPM or >120 FPM</i> | <i>Red LED and continuous audible alarm</i> |
| | Velocity Low | <i><90 FPM</i> | <i>Yellow LED and pulsing audible alarm</i> |
| | Velocity High | <i>>110 FPM</i> | <i>Yellow LED and pulsing audible alarm</i> |

face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pushing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm.

FOUR BENCH FUME HOODS & DUAL SUPPLY SINGLE GENERAL EXHAUST PRESSURE/TEMPERATURE CONTROL

| PROJECT: SAMPLE Standard Application #27 | | | | | | | | | | PROJECT NO.: | | | | DATE |
|------------------------------------------|--------------------|------------------------|----------------|---------------|---------|---------|---------------|------------------|-------------|--------------|--|--|--|------|
| ROOM NUMBER | DESCRIPTION | FUNCTION | PART NUMBER | TAG # | MAX CFM | MIN CFM | FAIL POSITION | RANGE | COMMENTS | | | | | |
| APP27 | Controller | Room Control | LC-2 | LC2-APP27 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Damper & Actuator | Room Supply #1 Control | CPP-P-14-NC-IL | CPPS-APP27-S1 | 2250 | 220 | N/C | 4-20 mA | | | | | | |
| APP27 | DP Transmitter | Room Supply #1 CFM | FTR-0.5 | FTRS-APP27-S1 | 3030 | 0 | n/a | 0.0 to 0.5"wc | | | | | | |
| APP27 | Damper & Actuator | Room Supply #2 Control | CPP-P-14-NC-IL | CPPS-APP27-S2 | 2250 | 220 | N/C | 4-20 mA | | | | | | |
| APP27 | DP Transmitter | Room Supply #2 CFM | FTR-0.5 | FTRS-APP27-S2 | 3030 | 0 | n/a | 0.0 to 0.5"wc | | | | | | |
| APP27 | Transmitter | Room Pressure | SPR-1A | SPR-APP27 | n/a | n/a | n/a | 0.0 to 0.01"wc | | | | | | |
| APP27 | Press. Sensor Pair | Room Pressure | SPS-2T | SPS-APP27-1A | n/a | n/a | n/a | | | | | | | |
| | | | | SPS-APP27-1B | n/a | n/a | n/a | | | | | | | |
| APP27 | Temp. Sensor | Space Temperature | STAT | STAT-APP27 | n/a | n/a | n/a | -32 to 132 *F | | | | | | |
| APP27 | Damper & Actuator | General Exhaust | PDA-10-NO-IL | PDAX-APP27 | 1090 | 0 | N/O | 4-20 mA | | | | | | |
| APP27 | Controller | Hood #1 Control | HC1 | HC1-APP27-1 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Transmitter | Hood #1 Face Velocity | FVR-1A | FVR-APP27-1 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP27 | Air Foil Pitot | Hood #1 Face Velocity | AFP2-60 | AFP-APP27-1 | n/a | n/a | n/a | n/a | 6-foot Hood | | | | | |
| APP27 | Tubing Kit | Hood #1 Face Velocity | TBK-1 | TBK-APP27-1 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Damper & Actuator | Hood #1 Exhaust | PDA-10-NO-HM | PDAX-APP27-1 | 1250 | 250 | N/O | 4-20 mA | 6-foot Hood | | | | | |
| APP27 | Controller | Hood #2 Control | HC1 | HC1-APP27-2 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Transmitter | Hood #2 Face Velocity | FVR-1A | FVR-APP27-2 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP27 | Air Foil Pitot | Hood #2 Face Velocity | AFP2-48 | AFP-APP27-2 | n/a | n/a | n/a | n/a | 5-foot Hood | | | | | |
| APP27 | Tubing Kit | Hood #2 Face Velocity | TBK-1 | TBK-APP27-2 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Damper & Actuator | Hood #2 Exhaust | PDA-8-NO-HM | PDAX-APP27-2 | 1000 | 200 | N/O | 4-20 mA | 5-foot Hood | | | | | |
| APP27 | Controller | Hood #3 Control | HC1 | HC1-APP27-3 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Transmitter | Hood #3 Face Velocity | FVR-1A | FVR-APP27-3 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP27 | Air Foil Pitot | Hood #3 Face Velocity | AFP2-60 | AFP-APP27-3 | n/a | n/a | n/a | n/a | 6-foot Hood | | | | | |
| APP27 | Tubing Kit | Hood #3 Face Velocity | TBK-1 | TBK-APP27-3 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Damper & Actuator | Hood #3 Exhaust | PDA-10-NO-HM | PDAX-APP27-3 | 1250 | 250 | N/O | 4-20 mA | 6-foot Hood | | | | | |
| APP27 | Controller | Hood #4 Control | HC1 | HC1-APP27-4 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Transmitter | Hood #4 Face Velocity | FVR-1A | FVR-APP27-4 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP27 | Air Foil Pitot | Hood #4 Face Velocity | AFP2-48 | AFP-APP27-4 | n/a | n/a | n/a | n/a | 5-foot Hood | | | | | |
| APP27 | Tubing Kit | Hood #4 Face Velocity | TBK-1 | TBK-APP27-4 | n/a | n/a | n/a | n/a | | | | | | |
| APP27 | Damper & Actuator | Hood #4 Exhaust | PDA-8-NO-HM | PDAX-APP27-4 | 1000 | 200 | N/O | 4-20 mA | 5-foot Hood | | | | | |
| | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | |

NOTE:
SEE DRAWINGS A950051, A9500511 &
A9500515 FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.


The diagram illustrates a fume exhaust system for four 'AUTO-FLUX' units. The system includes a supply line with a pressure-reducing valve (PRV) and a check valve (C), a main supply line with a pressure-reducing valve (PRV) and a check valve (C), and a fume exhaust line. Each unit has a fume exhaust port labeled 'EXH' and a water trap labeled 'W/T'. The units are labeled 'AUTO-FLUX'.

| | | |
|----------------------------------------|-----------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| TYPICAL FUME HOOD LAB APPLICATION H 27 | DRAWN LIP ENGR. REV. CHAS. RUS. DATE: 9-17-92 SHEET | AUTO-FLOW PA 1086 BY S.S. MODEL PA 1432 MAX. APPLICABLE CONTAMINANT NOT TO EXCEED: PRESET FLOW CONTROL: 1.00 MGD DRAWING NO. A9400427 |
|----------------------------------------|-----------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|

AUTO-ELDW
ONE HONOLULU DRIVE, HONOLULU, HI 96813
PHONE: 808/591-1111
FAX: 808/591-1112
E-MAIL: info@auto-eldw.com
WWW: www.auto-eldw.com

Application #28

Dual Supply Single General
Exhaust for Air Flow Control
with Full Temperature Control
and Four Bench Hoods



In this Sequence:

LABORATORY

- *General*
- *Air Flow Control*
- *Temperature Control*
- *Volume Limit Control*
- *Emergency Mode*
- *Alarming*

FUME HOODS

- *General*
- *Control*
- *Exhaust Air Flow
Monitoring*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #28: Dual Supply Single General Exhaust for Air Flow Control with Full Temperature Control and Four Bench Hoods

LABORATORY

General

The laboratory is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on laboratory pressure control.

Air Flow Control

The controller measures the laboratory supply air flow with a parallel plate pitot array in both of the laboratory supply air ducts and measures the exhaust air flow leaving the space with a parallel plate pitot array in the laboratory general exhaust duct. These measurements are scaled into cubic feet per minute (CFM). The controller executes its process loops using proportional, integral, and derivative control action. The supply air flows are the measured variables for two of the controller's PID loops. The setpoint for each supply air flow is half of the minimum air flow rate required for the desired air changes per hour, reset from the temperature control signal or half of the air flow rate required to satisfy the fume hood makeup air volume, whichever is greater. The supply makeup air flow rate is determined by the following equation:

Supply Makeup Flow Setpoint CFM = (Total Hood Exhaust Air Flow CFM + Minimum General Exhaust Air Flow CFM - Desired Differential CFM Setpoint)

The controller's PID outputs are used to control output to each of the laboratory supply air dampers. The laboratory supply air dampers are positioned as required by the controller with electronic damper

actuators. The damper actuators incorporate internal feedback for positive positioning. The controller executes its control loop five times per second. When supply air flow falls below the desired setpoint the controller modulates the normally closed laboratory supply air dampers toward the open position. When supply air flow rises above the desired setpoint the controller modulates the normally closed laboratory supply air dampers toward the closed position.

The exhaust air flow setpoint is based on a calculation derived to maintain a constant differential CFM. This differential CFM is the mathematical result of the total exhaust air leaving the space minus the total supply air entering the space. The differential CFM setpoint is a constant value that is user adjustable. The exhaust air flow setpoint is based on the following equation:

Exhaust CFM Setpoint = Current Supply CFM + Δ CFM Setpoint - Total Hood Exhaust CFM

The exhaust air flow PID is controlled to maintain this differential CFM. The controller's PID output is used to control output to the laboratory general exhaust air damper. The laboratory general exhaust air damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop five times per second. When general exhaust air flow falls below the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the open position. When general exhaust air flow rises above the desired setpoint the controller modulates the normally open laboratory general exhaust air damper toward the closed position.

Temperature Control

The controller measures the temperature of the supply air in the supply duct serving the laboratory space and the temperature of the laboratory space itself. These temperatures combined with the measured supply air flow rate are used to calculate the instantaneous BTUs used by the space to maintain its temperature. The equation for this process is as follows:

Sequence of Operation

$BTUs = [(Supply\ Discharge\ Temperature - Space\ Temperature) \div Supply\ Air\ Flow\ Rate] \div 1.085$

This equation is calculated twice, one for each supply duct. One of the controller's PID loops is used in conjunction with the measured space temperature and its required setpoint to determine the desired instantaneous BTU setpoint. Half of this BTU setpoint is used in conjunction with the respective calculated current BTUs and a PID loop to control each of the laboratory supply duct reheat coils. Additionally, the BTU setpoint and calculated current BTUs work in a reset fashion with the minimum supply flow setpoint and each of the associated PID loops controlling the laboratory supply air dampers. If the BTUs fall below the required setpoint the controller decreases the minimum supply air flow setpoint associated with the respective PID loop controlling each of the supply air dampers toward its minimum value (the minimum design CFM for the required air change rate). As the calculated BTUs continue to fall below the required setpoint the controller modulates the outputs to the reheat coil valves toward the open position. If calculated BTUs rise above the desired setpoint the controller modulates the reheat coil outputs toward the closed position. On a continued rise above the desired setpoint the controller resets the minimum supply air flow setpoint associated with the respective PID loop controlling each of the supply air dampers toward its maximum value (the maximum design CFM for the lab supply).

Volume Limit Control

Should the measured supply air flow volume fall below its desired minimum air flow rate (adjustable) as required to maintain the necessary air changes per hour, each of the supply flow control PID loops is overridden and the laboratory supply dampers are controlled to maintain the desired minimum air flow. Should the measured supply air flow volume rise above its desired maximum air flow rate (adjustable) each of the supply flow control PID loops is overridden and the laboratory supply dampers are controlled to maintain the desired maximum air flow.

Emergency Mode

The controller has the capability to accept an external contact from an emergency push button which indicates the emergency purge mode. In this mode the laboratory general exhaust air and supply

air dampers are overridden to a fixed (user-adjustable) position and an alarm is transmitted over the communications network to the remote operator interface device. These dampers remain in this position until the emergency push button is reset to its normal state.

Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power or actuator position feedback not meeting its commanded position. Additionally, the controller also allows for alarms on the supply air flows, general exhaust air flow, supply duct temperatures and space temperature measurements based on user-adjustable limits. These alarms are not annunciated locally. They are transmitted over the communications network to the remote operator interface device.

FUME HOODS

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has built-in local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LED(s) capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading acts as the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, which is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|----------------|------------------------|----------------------------------|---------------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm. When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If the face velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for use by the laboratory pressure controller in calculating total lab exhaust volume. The controller provides a dedicated input and output for receiving a cascaded flow signal from another fume hood controller, adding its own exhaust flow value and transmitting a totalized signal to either another fume hood controller or the laboratory pressure controller.

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.

FOUR BENCH FUME HOODS & DUAL SUPPLY SINGLE GENERAL EXHAUST AIR FLOW/TEMPERATURE CONTROL

| PROJECT: SAMPLE Standard Application #28 | | | | | | | | | | PROJECT NO. : | | | | DATE |
|------------------------------------------|-------------------|-----------------------|----------------|---------------|---------|---------|---------------|------------------|--|---------------|--|--|--|-------------|
| ROOM NUMBER | DESCRIPTION | FUNCTION | PART NUMBER | TAG # | MAX CFM | MIN CFM | FAIL POSITION | RANGE | | | | | | COMMENTS |
| APP28 | Controller | Room Control | LC-2 | LC2-APP28 | n/a | n/a | n/a | n/a | | | | | | 05/06/2000 |
| APP28 | Damper & Actuator | Room Supply #1 CFM | CPP-P-14-NC-IL | CPPS-APP28-S1 | 2000 | 250 | N/C | 4-20 mA | | | | | | |
| APP28 | DP Transmitter | Room Supply #1 CFM | FTR-0.25 | FTRS-APP28-S1 | 2140 | 0 | n/a | 0.0 to 0.25"wc | | | | | | |
| APP28 | Damper & Actuator | Room Supply #2 CFM | CPP-P-14-NC-IL | CPPS-APP28-S2 | 2000 | 250 | N/C | 4-20 mA | | | | | | |
| APP28 | DP Transmitter | Room Supply #2 CFM | FTR-0.25 | FTRS-APP28-S2 | 2140 | 0 | n/a | 0.0 to 0.25"wc | | | | | | |
| APP28 | Damper & Actuator | Room Gen. Exh. CFM | CPP-P-10-NO-IL | CPPX-APP28 | 1090 | 0 | N/C | 4-20 mA | | | | | | |
| APP28 | DP Transmitter | Room Gen. Exh. CFM | FTR-0.25 | FTRX-APP28 | 1092 | 0 | n/a | 0.0 to 0.25"wc | | | | | | |
| APP28 | Temp. Sensor | Space Temperature | STAT | STAT-APP28 | n/a | n/a | n/a | -32 to 132 °F | | | | | | |
| APP28 | Controller | Hood #1 Control | HC1 | HC1-APP28-1 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | DP Transmitter | Hood #1 Face Velocity | FVR-1A | FVR-APP28-1 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP28 | Air Foil Pitot | Hood #1 Face Velocity | AFP2-36 | AFP-APP28-1 | n/a | n/a | n/a | n/a | | | | | | 4-foot Hood |
| APP28 | Tubing Kit | Hood #1 Face Velocity | TBK-1 | TBK-APP28-1 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | Damper & Actuator | Hood #1 Exhaust | CPP-P-8-NO-HM | CPPX-APP28-1 | 750 | 150 | N/O | 4-20 mA | | | | | | 4-foot Hood |
| APP28 | DP Transmitter | Hood #1 Exhaust CFM | FTR-0.5 | FTRX-APP28-1 | 988 | 0 | n/a | 0.0 to 0.5"wc | | | | | | |
| APP28 | Controller | Hood #2 Control | HC1 | HC1-APP28-2 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | DP Transmitter | Hood #2 Face Velocity | FVR-1A | FVR-APP28-2 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP28 | Air Foil Pitot | Hood #2 Face Velocity | AFP2-60 | AFP-APP28-2 | n/a | n/a | n/a | n/a | | | | | | 6-foot Hood |
| APP28 | Tubing Kit | Hood #2 Face Velocity | TBK-1 | TBK-APP28-2 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | Damper & Actuator | Hood #2 Exhaust | CPP-P-10-NO-HM | CPPX-APP28-2 | 1250 | 250 | N/O | 4-20 mA | | | | | | 6-foot Hood |
| APP28 | DP Transmitter | Hood #2 Exhaust CFM | FTR-0.5 | FTRX-APP28-2 | 1544 | 0 | n/a | 0.0 to 0.5"wc | | | | | | |
| APP28 | Controller | Hood #3 Control | HC1 | HC3-APP28-3 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | DP Transmitter | Hood #3 Face Velocity | FVR-1A | FVR-APP28-3 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP28 | Air Foil Pitot | Hood #3 Face Velocity | AFP2-36 | AFP-APP28-3 | n/a | n/a | n/a | n/a | | | | | | 4-foot Hood |
| APP28 | Tubing Kit | Hood #3 Face Velocity | TBK-1 | TBK-APP28-3 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | Damper & Actuator | Hood #3 Exhaust | CPP-P-8-NO-HM | CPPX-APP28-3 | 750 | 150 | N/O | 4-20 mA | | | | | | 4-foot Hood |
| APP28 | DP Transmitter | Hood #3 Exhaust CFM | FTR-0.5 | FTRX-APP28-3 | 988 | 0 | n/a | 0.0 to 0.5"wc | | | | | | |
| APP28 | Controller | Hood #4 Control | HC1 | HC1-APP28-4 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | DP Transmitter | Hood #4 Face Velocity | FVR-1A | FVR-APP28-4 | n/a | n/a | n/a | 0.0 - 0.0015 "wc | | | | | | |
| APP28 | Air Foil Pitot | Hood #4 Face Velocity | AFP2-60 | AFP-APP28-4 | n/a | n/a | n/a | n/a | | | | | | 6-foot Hood |
| APP28 | Tubing Kit | Hood #4 Face Velocity | TBK-1 | TBK-APP28-4 | n/a | n/a | n/a | n/a | | | | | | |
| APP28 | Damper & Actuator | Hood #4 Exhaust | CPP-P-10-NO-HM | CPPX-APP28-4 | 1250 | 250 | N/O | 4-20 mA | | | | | | 6-foot Hood |
| APP28 | DP Transmitter | Hood #4 Exhaust CFM | FTR-0.5 | FTRX-APP28-4 | 1544 | 0 | n/a | 0.0 to 0.5"wc | | | | | | |
| | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | |

The diagram illustrates a four-burner gas cooktop with a built-in fume exhaust system. Each burner is equipped with an 'AUTO-FLAME' control knob. A central fume exhaust duct is connected to the cooktop, leading to a vertical duct and then to a horizontal duct that exits through a wall. The diagram includes labels for 'FUME EXHAUST' and 'EXH'.

FUME EXHAUST


STAT

TYPICAL FUME HOOD LAB
APPLICATION # 28

4 BENCH TYPE FUME HOODS

Application #29

Constant Volume On/Off or
Switched Setpoint Control



In this Sequence:

- *General*
- *Air Flow Control*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #29: Constant Volume On/Off or Switched Setpoint Control

General

The constant volume device is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions independently of the network for all lab control aspects so that a failure of the communications network has no effect over control.


Air Flow Control

The controller measures the flow through the device with a parallel plate pitot array in the duct servicing the device. These measurements are scaled into cubic feet per minute. The controller executes its process loops using proportional, integral, and derivative control action. The controller and sensor combination act as a completely pressure independent control unit. Air flow is the measured variable for the controller's PID loop. The setpoint for air flow is fixed and user adjustable. The setpoint has the capability to change to a different fixed (user-adjustable) value based on input to the controller from an external switch or based on a time of day schedule. The controller's PID output is used to control output to the device's air regulating damper. The air regulating damper is positioned as required by the controller with an electronic damper actuator. The damper actuators incorporate internal feedback for positive positioning. The controller executes its control loop of five times per second. If air flow falls below the desired setpoint the controller modulates the air regulating damper toward the open position. If the air flow rises above the desired setpoint the controller modulates the air regulating damper toward the closed position.

Alarming

The controller provides the capability to alarm based on an actuator failure related to either loss of power

or actuator position feedback that does not meet the commanded position. The controller also allows for alarms on the air flow measurements based on user-adjustable limits. These alarms annunciate locally or remotely or a combination of both.




BENCH TYPE FUME HOOD
W/EXHAUST FLOW MEASUREMENT

| | | | | | | | |
|----------------------------|--|------------------|--|---------------|--|----------------------|--|
| DRAWN: LUP | | CHECKED: RVT | | DATE: 9-13-95 | | SHEET: 1 | |
| TYPICAL FUME HOOD LAB | | APPLICATION # 29 | | AUTOFLOW | | DRAWING NO. A9400429 | |
| BY: J. J. JONES, P.E. | | DATE: 9-13-95 | | SHEET: 1 | | DRAWING NO. A9400429 | |
| CONTAINS: FUME HOOD | | DATE: 9-13-95 | | SHEET: 1 | | DRAWING NO. A9400429 | |
| CONTAINS: VOLUME OF FUMES | | DATE: 9-13-95 | | SHEET: 1 | | DRAWING NO. A9400429 | |
| CONTAINS: SETPOINT CONTROL | | DATE: 9-13-95 | | SHEET: 1 | | DRAWING NO. A9400429 | |

Application #30

Glove Box Pressure Control



In this Sequence:

- *General*
- *Pressure Control*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW®

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #30: Glove Box Pressure Control

General


The glove box is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller is capable of interfacing with local network interface devices such as hand-held or wall-mounted local displays. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on control.

Pressure Control

The total pressure differential is measured between the glove box interior and the surrounding space. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is $-0.01''\text{wc}$ with respect to the surrounding space, a setpoint that is adjustable either locally or remotely. The controller's PID output is used to control output to the glove box exhaust air damper. The exhaust air damper is positioned by the controller with an electronic damper actuator. The damper actuators incorporate internal feedback for positive positioning. The controller executes its control loop five times per second. If the pressure differential falls below the desired setpoint the controller modulates the exhaust air damper toward the open position. If the pressure differential goes above the desired setpoint the controller modulates the exhaust air damper toward the closed position.

Alarming

The controller provides the capability to alarm based on actuator failures related to either loss of power or actuator position feedback that does not meet commanded position. The controller also allows for alarms of the differential pressure measurements based on user-adjustable limits. These alarms announce both locally and remotely.




NOTE:
SEE DRAWING A950057
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | | | | | |
|-----------------------|--|----------------|--|----------------------|--|----------|--|
| DRAWN: LIP | | CHECKED: RWT | | DATE: 5-12-95 | | SHEET: 1 | |
| TYPICAL FUME HOOD LAB | | APPROVED: RWT | | DATE: 5-12-95 | | SHEET: 1 | |
| APPLICATION #30 | | PROJECT: CHDRA | | DRAWING NO: A9400430 | | AUTO-FLD | |

Application #31

Biological Safety Cabinet



In this Sequence:

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW®

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #31: Single Bench Biological Safety Cabinet

General

The biological safety cabinet is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LEDs capable of displaying current cabinet operating status in up to three different conditions
- audible alarm
- audible alarm silence
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on biological safety cabinet control.

Control

The controller measures the exhaust air flow volume leaving the biological safety cabinet with a parallel plate array in the cabinet exhaust duct. This measurement is scaled into cubic feet per minute (CFM) for use by the biological safety cabinet controller. The controller executes a process loop using proportional, integral, and derivative control action. The exhaust air flow reading is the measured variable for the controller's PID loop. The setpoint for exhaust air flow is a flow rate that translates to an average face velocity of 75 (Class I, Type A) or 100 (Class I, Type B1, 2, &3) feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software

providing an equal percentage control output to the biological safety cabinet exhaust damper. The biological safety cabinet exhaust flow damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When exhaust flow falls below the desired setpoint the controller modulates the normally open biological safety cabinet exhaust damper toward the open position. If exhaust flow continues to fall below setpoint, the controller generates a local alarm (as described in the alarming section) after an adjustable time delay. When exhaust flow rises above the desired setpoint the controller modulates the normally open biological safety cabinet exhaust damper toward the closed position. If exhaust flow continues to rise above setpoint, the controller generates a local alarm after an adjustable time delay.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the cabinet. Initiation of either device positions the biological safety cabinet exhaust damper to an user-defined fixed position. In this mode the controller display indicates that the hood is in emergency mode while the audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally or a reset of the remote device.


Alarming

The controller has the ability to initiate alarms both locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | |
|------------------------|--------------------|--------------------------------------|
| Velocity Unsafe | <20 FPM or >20 FPM | Red LED and continuous audible alarm |
| Velocity Low | <10 FPM | Yellow LED and pulsing audible alarm |
| Velocity High | >10 FPM | Yellow LED and pulsing audible alarm |

Table 1

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms can be silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that mimics the operation of the built-in audible alarm.



NOTE:
SEE DRAWING A950058
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.


| | | | |
|---------------|--|--------------------------------|--|
| DRAWN BY | | AUTO-FLY | |
| ENGR. BY | | ON 10/10/87 BY TUPPEL, PA 1932 | |
| CHECKED BY | | BIO SAFETY HOOD | |
| DATE: 8-13-95 | | FLOW DIAGRAM | |
| SHEET | | DRAWING NO. A9400431 | |

TYPICAL FUME HOOD LAB
APPLICATION H 31

BIO SAFETY HOOD

Application #32

Single Bench Fume Hood with
Exhaust Bypass Control



In this Sequence:

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #32: Single Bench Fume Hood with Exhaust Bypass Control

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LEDs capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Fume Hood Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an equal percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If velocity continues to fall below setpoint, the controller

generates a local alarm (as described in the alarming section) after an adjustable time delay. When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If velocity continues to rise above setpoint, the controller generates a local alarm after an adjustable time delay.

Exhaust Bypass Damper Control

The controller measures static pressure in the common exhaust duct plenum using a duct static pressure sensor. The controller executes a process loop using proportional, integral, and derivative control action. The exhaust duct static pressure reading is the measured variable for the controller's PID loop. The setpoint for exhaust duct static pressure is determined in the field and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an equal percentage control output to the exhaust bypass damper. The exhaust bypass damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. When exhaust duct static pressure falls below the desired setpoint the controller modulates the normally closed exhaust bypass exhaust damper toward the open position.. When exhaust duct static pressure rises above the desired setpoint the controller modulates the normally closed exhaust bypass damper toward the closed position.

Emergency Override


The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood.. Initiation of either device positions the hood exhaust damper and exhaust bypass damper to a user-defined fixed position. In this mode the controller display indicates that the hood is in emergency mode while the audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally.

Alarming

The controller has the ability to initiate alarms both locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1) :

| | | | |
|----------------|------------------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms can be silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that mimics the operation of the built-in audible alarm.




NOTE:
SEE DRAWINGS A950051 & A950051S
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

BENCH TYPE FUME HOOD

| | | | |
|----------------------------------------|---------------|-----------------------------------|---------------|
| TYPICAL FUME HOOD LAB APPLICATION # 32 | | AUTO-FLUX | |
| DESIGN: LIP | ENGR: RVT | REVISION: 1 | DATE: 10/1/82 |
| CHECK: RVT | DATE: 1/13/93 | SINGLE HOOD - BYPASS FLOW DIAGRAM | |
| SHEET 1 | | DRAWING NO. A9400432 | |

Application #33

Single Bench Fume Hood with
General Exhaust Air Flow
Tracking



In this Sequence:

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #33: Single Bench Fume Hood with General Exhaust Air Flow Tracking

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LEDs capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Fume Hood Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an adjusted percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller using an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If flow continues to fall below setpoint, after an adjustable time delay the controller generates a local alarm (as described in the

alarming section). When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If velocity continues to rise above setpoint the controller generates a local alarm after an adjustable time delay.

Fume Hood Exhaust Air Flow Monitoring

The controller measures the exhaust air flow volume leaving the fume hood using a parallel plate pitot array in the fume hood exhaust duct. This measurement is scaled into cubic feet per minute for use by the controller in calculating total lab exhaust volume.

General Exhaust Air Flow Control

The controller measures the exhaust air flow volume leaving the general exhaust grille with a parallel plate pitot array and differential pressure transmitter in the general exhaust duct. This measurement is scaled into cubic feet per minute. The controller executes its process loop with proportional, integral, and derivative control action. The general exhaust air flow is the measured variable for the controller's PID loop. The general exhaust air flow setpoint is based on the following equation:

General Exhaust Setpoint = (Total exhaust flow capacity- fume hood exhaust flow)

The setpoint is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory general exhaust air damper. The general exhaust air damper is positioned as required by the controller through an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. If general exhaust air flow falls below the desired setpoint the controller modulates the normally closed general exhaust air damper toward the open position. On a rise in general exhaust air flow above the desired setpoint the controller modulates the normally closed general exhaust air damper toward the closed position.

Sequence of Operation



Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper to a user-defined fixed position. In emergency mode the controller display audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally, or a resetting of the remote device.


After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that parallels the operation of the built-in audible alarm. The controller also provides an alarm for the exhaust flow value based on user adjustable limits. This alarm does not annunciate locally, it is transmitted over the communications network to the remote operator interface device.

Alarming

The controller has the ability to initiate alarms locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1):

| | | | |
|---------|-----------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms are silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period.




NOTE:
SEE DRAWINGS A950053 & A9500516
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

BENCH TYPE FUME HOOD
W/EXHAUST FLOW MEASUREMENT

| | | |
|----------------------------------------|-----------|------------------|
| DESIGN LIP | AUTO-FLOW | EXHAUST, PA 1532 |
| ENGR. BY | EXHAUST | EXHAUST |
| CHKD. BY | EXHAUST | EXHAUST |
| DATE: 11/17/95 | EXHAUST | EXHAUST |
| SHEET | EXHAUST | EXHAUST |
| TYPICAL FUME HOOD LAB APPLICATION # 33 | | |
| DRAWING NO. A9400433 | | |

Application #34

Isolation Room Control



In this Sequence:

- *General*
- *Configurationl*
- *Control*
- *Alarming*

Bill of Material
Drawings

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #34: Isolation Room Control

General

The isolation room is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities so that it can interface with hand-held or wall-mounted local display devices. The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all lab control aspects so that a failure of the communications network has no effect on fume hood control.

Configuration

The isolation room ventilation is configured as described in this section. The supply air to the space is controlled separately to its design requirements. The space is served by an exhaust air duct with a modulating control damper. A space pressure sensing device measures the total pressure differential across the barrier wall between the controlled space and the adjacent space. The controller panel includes a face-mounted alarm light, alarm horn, horn silence button, space pressure LCD display, isolation/normal selector switch, and a positive/negative selector switch.

Mode Control

When the isolation/normal selector switch is in the normal position the exhaust control damper is modulated to maintain the space pressure at the neutral setpoint. When the isolation/normal selector switch is in the isolation position the exhaust control damper is modulated to maintain the space pressure at the setpoint determined by the position of the positive/negative selector switch. When the positive/negative selector switch is in the positive position the space pressure setpoint is a user-adjustable positive value and the space pressure sensor is set to measure a positive differential in the space with reference to the adjacent area. When the positive/negative selector switch is in the negative position the space pressure setpoint is a user-adjustable negative value and the space pressure



sensor is set to measure a negative differential in the space with reference to the adjacent area.

Pressure Control

The total pressure differential is measured between the isolation room and the adjacent space. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The controller's PID output is used to control either the exhaust control damper or the return air control damper based on the current operating mode. The control damper is positioned as required by the controller with a fully electronic damper actuator. The controller executes its complete control loop a minimum of five times per second.

Alarming


The controller provides the capability to detect an alarm based on an actuator failure related to either loss of power or actuator position feedback that does not meet commanded position. The controller also allows for alarming of the differential pressure measurements based on user-adjustable limits. The alarm panel light and horn energize whenever the space pressure is outside the user-adjustable limits for the current setpoint. When a pressure alarm occurs the user has the option of silencing the audible alarm horn by pressing the alarm silence button on the panel face. The pressure display continuously shows the current space pressure value.



NOTE:
SEE DRAWINGS A9500515 & A9500521
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|------------------|----------------------|---------------------------------|--|
| ISOLATION ROOM | | AUTO-FLW | |
| APPLICATION # 34 | | ONE THERMOSTAT PER ROOM, PA IN2 | |
| DESIGN: R/V | ENGR: R/V | SINGLE SUPPLY/EXHAUST | |
| CHKD: R/V | DATE: 9-11-95 | PRESSURE FLOW DIAGRAM | |
| SHEET | DRAWING NO. A9400434 | | |

■ (MOUNT PANEL WITHIN 12' OF SPACE PRESSURE SENSORS)




FIELD PANEL # 1
FP-1

| | | | |
|-----------------------|-----------|-----------------------|--------------|
| ISOLATION ROOM | | AUTO-FLY | |
| APPLICATION # 34 | | DRAWING NO. A9400434A | |
| DESIGN LIP | ENGR. RVT | DATE: 11-17-75 | BY: J. J. J. |
| CHG. RVT | CHG. RVT | DATE: 11-17-75 | BY: J. J. J. |
| SHEET SUPPLEMENT | | SHEET SUPPLEMENT | |
| PRESSURE FLOW DIAGRAM | | PRESSURE FLOW DIAGRAM | |

Application #35

Single Bench Fume Hood with
Single Supply Pressure
Control



In this Sequence:

- *General*
- *Control*
- *Emergency Override*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #35: Single Bench Fume Hood with Single Supply Pressure Control

General

The fume hood is equipped with a stand-alone networkable microprocessor-based controller. The controller has password protection at both the local and remote interface levels. This protection supports multiple levels of access. The controller has local operator interface capabilities consisting of but not limited to the following:

- two line LCD or LED display capable of simultaneously displaying scaled point values and English language text point descriptors
- LEDs capable of displaying current hood operating status in up to three different conditions
- audible alarm
- audible alarm silence key
- key style interface to enter password and gain menu access

The networking capabilities allow for remote operator interface as well as data/alarm logging and storage. The controller functions as a stand-alone unit for all hood control aspects so that a failure of the communications network has no effect on fume hood control.

Fume Hood Control

The controller measures face velocity directly in the plane of the fume hood sash opening using the air foil pitot method. The controller executes a process loop using proportional, integral, and derivative control action. The face velocity reading is the measured variable for the controller's PID loop. The setpoint for face velocity is usually set as 100 feet per minute, and is adjustable either locally or remotely. The controller's PID output is calculated by the controller software providing an equal percentage control output to the fume hood exhaust damper. The fume hood exhaust damper is positioned as required by the controller with an electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its control loop twenty times per second. When face velocity falls below the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the open position. If velocity continues to fall below setpoint, the controller

generates a local alarm (as described in the alarming section) after an adjustable time delay. When face velocity rises above the desired setpoint the controller modulates the normally open fume hood exhaust damper toward the closed position. If velocity continues to rise above setpoint, the controller generates a local alarm after an adjustable time delay, as described in the alarming section.

Lab Pressure Control

The total pressure differential is measured between the laboratory space and the adjacent corridor. The controller executes a process loop using proportional, integral, and derivative control action. The total pressure differential reading is the measured variable for the controller's PID loop. The setpoint for total pressure differential is - 0.005"wc with respect to the adjacent corridor, and is adjustable either locally or remotely. The controller's PID output is used to control output to the laboratory supply damper. The laboratory supply damper is positioned by the controller with a fully electronic damper actuator. The damper actuator incorporates internal feedback for positive positioning. The controller executes its complete control loop five times per second. If the total pressure differential falls below the desired setpoint the controller modulates the normally closed laboratory supply damper toward the open position. If total pressure differential goes above the desired setpoint the controller modulates the normally closed laboratory supply damper toward the closed position.

Emergency Override

The controller has an emergency override key on its local operator interface. In addition, the controller has an external emergency override input to cover initiating devices located away from the hood. Initiation of either device positions the hood exhaust damper and the lab supply damper to user-defined fixed positions. In this mode the controller display indicates that the hood is in emergency mode while the audible alarm sounds and the red LED pulses. Clearing this mode requires that the controller display's enter key be pushed if the condition is initiated locally or reset of the remote device.

Sequence of Operation




Fume Hood Alarming


The controller has the ability to initiate alarms both locally and remotely (if required) based on the following (adjustable) face velocity values (See Table 1) :

| | | | |
|---------|-----------------|---------------------|--------------------------------------|
| Table 1 | Velocity Unsafe | <80 FPM or >120 FPM | Red LED and continuous audible alarm |
| | Velocity Low | <90 FPM | Yellow LED and pulsing audible alarm |
| | Velocity High | >110 FPM | Yellow LED and pulsing audible alarm |

In addition, the controller provides the capability to alarm based on an actuator failure related either to loss of power or because actuator position feedback did not meet commanded position. In this alarm condition the controller display's red LED and the audible alarm pulse. Audible alarms can be silenced by pressing the silence key on the controller's display module. This key does not cancel the alarm condition but does silence the alarm for a programmable time period. After the time delay the audible alarm returns if the condition has not been corrected. The controller provides a dedicated digital output for initiating a remote alarm klaxon or remote alarm light that mimics the operation of the built-in audible alarm.



SUPPLY DAMPER CONTROLLED BY HC1 CONTROLLER.




NOTE:
SEE DRAWINGS A950051 & A9500520
FOR SPECIFIC WIRING DETAILS
AND PART NUMBER INFORMATION.

| | | | |
|---------------------------------|-----------|-------------------------------|-------|
| TYPICAL FUME HOOD LAB | | DRAWING NO. A9400435 | |
| APPLICATION # 35 | | DRAWING NO. A9400435 | |
| DESIGN LIP | ENGINEER | DATE | SHEET |
| ENGR. BY | CHIEF. BY | DATE: 8-12-95 | 1 |
| SINGLE SUPPLY SINGLE BENCH HOOD | | PRESSURE CONTROL FLOW DIAGRAM | |
| AUTO-FLOW | | MODEL, PA 1942 | |

BENCH TYPE FUME HOOD

Application #36

Single Bench Fume Hood-
Face Velocity Monitoring
Only



In this Sequence:

- *General*
- *Hood Monitoring*
- *Alarming*

Bill of Material
Drawing

AUTO-FLOW[®]

Sequence of Operation

Laboratory and Fume Hood Controls Engineering Guide

Application #36: Single Bench Fume Hood- Face Velocity Monitoring Only

General

The fume hood is equipped with an electronic fume hood monitor. The monitor has built-in local operator interface capabilities that include the following:

- Normal flow indicator
- Alarm flow indicator
- Test/reset button
- Power on indicator
- Audible alarm
- Alarm setpoint adjustment potentiometer


Hood Monitoring

The monitor measures face velocity by the flow of air passing through the air inlet located on the face of the HM1. The air flow passes over dual hot/cold thermistors and exits through the sidewall of the hood.

Alarming

The controller has the ability to initiate alarms both locally and remotely (if required) based on the following (adjustable) face velocity values:

| | | |
|----------------------|-----------------------------------------|----------------------------------------------------------------------------------------------|
| Velocity Safe | >80 FPM (adjustable from 50 to 150 FPM) | Green LED on |
| Velocity Low | <80 FPM (adjustable from 50 to 150 FPM) | Red LED and continuous audible alarm <i>OPTIONAL</i> NC relay contact (24V 0.3A DC or AC) |



BENCH TYPE FUME HOOD MONITOR ONLY

| | | | |
|-----------------------|---------------|-----------------------------|----------|
| TYPICAL FUME HOOD LAB | | AUTO-FLOW | |
| APPLICATION #36 | | DATE: 10/10/83 BY: J. J. J. | |
| DESIGNER: RWT | ENGINEER: RWT | DATE: 8-13-95 | SHEET: 1 |
| TYPICAL FUME HOOD | | FUME DIAGRAM | |
| DRAWING NO. A9400436 | | | |